

Electrophoresis Power Supply (QC-EPS-3000)



Main Features

1. Large-screen LCD displays voltage, current, power and timing time at the same time.
2. Voltage, current and power closed-loop control, realizing adjustment during operation.
3. With call recovery function.
4. After reaching the set time, it has the function of maintaining small current.
5. Perfect protection and early warning functions.
6. With memory storage function.
7. One machine with multiple slots, four parallel outputs.

Technical Specifications:

Voltage Output Range	2-300V (Resolution: 1V)
Current Output Range	5-2000mA (Resolution:1 mA)
Power Output Range	1-300W (Resolution:1W)
Input Voltage	220V, 50/60HZ
Parallel Output	4 groups
Timing Function	1min-99h59min, adjustable
External Size(L*W*H)	250*350*105 mm
Net/Gross Weight	3.5 kg/6 kg
Certificate	ISO & CE Certificate