

# Filtration System with Vacuum Pump



## Main Features

- One Year Warranty (QC 300)
- No air pollution, Maintenance free
- Clamp-free magnetic filter holder
- Drain out directly

### Oil-free & care-free:

- QC pumps are driven by piston that needs no lubricant, making them maintenance-free and oil-mist-free.

### Clamp-free with magnet:

- QC- 31 bottles are equipped with patented magnetic QC filter holders. Its strong magnetic seal enables one-hand and clamp-free operation.

### Drain out directly:

- QC- 31 can drain the filtrate away directly without disassembly.

### Technical Specifications:

Model	QC-FVP23
Max. power	65W
Max. current	0.7A
Max. vacuum	99 mbar
Max. flow rate	23 L/min
Motor rotation speed	1750 RPM
Horse power	1/8 HP
Noise level	50 dB
Hose barb	ID8 (5/16 inch)
Net weight	4.1 Kg
Dimension (LxWxH)	28.4 x 14.8 x 22.2 cm
Power	220 V
Vacuum Regulator	Yes
Overflow protection	Yes
Thermal protection	Yes
Certificate	ISO & CE Certificate