

Digital High-Speed Centrifuge Machine





Main Features

- Have the filtration and voltage regulation function.
- Run smoothly, precise separation.
- Small size, but large capacity.
- Various rotor, and easy to replace.
- High centrifugal power.
- Stop quickly.
- Super low working noise.
- Low vibration.

Technical Specification:

Model	QC-10Pro	QC-12Pro
Control mode	PWM electronic speed regulation	
Timing	1-99min/1-59sec	
Display	LCD	
Max relative centrifuge force	6780xg	9660xg
Max rotation speed	10000rpm±5%	12000rpm±5%
Max speed acceleration/reduction	≤12s / ≤15s	
Noise	≤47db(A)	
Input power	48W	
Shell protection level	AC100-240V / DC24V	
Allow surrounding temperature	5-40 °C / 80%RH	

Model	Rotor	Specifications
Multi Rotor		0.2/0.5/1.5/2ml micro tube*12 + 0.2ml 8 PCR stripe*4 + 0.5/1.0ml micro blood tube *12
Hematocrit Rotor		40mm capillary*12