

High Speed Mini Centrifuge (QC-MC-15Pro)



Main Features

1. Have the filtration and voltage regulation function.
2. Run smoothly, precise separation.
3. Small size, but large capacity.
4. Various rotor and easy to replace.
5. High centrifugal power.
6. Stop quickly.
7. Super low working noise.
8. Low vibration.

Technical Specifications:

Timing	10sec - 23h59min
Display	LCD Display
Max relative centrifuge force	14104xg
Max rotation speed/RPM	14500 RPM±5%
Power	55W
Shell protection level	AC100-240V / DC24V
Allow surrounding temperature	5-40°C / 80%RH
Dimensions	270*180*105 mm
Rotor	0.2/0.5/1.5/2.0ml centrifuge tube*12 0.5/1ml micro blood tube*12
Certificate	ISO & CE Certificate