

QC- Touch Screen Constant-Temperature Incubator



Main Features:

- Large LED touch screen.
- D type: with inner glass door.
- G type: with viewing window.
- Microcomputer controller.
- Stainless steel chamber.

Safety Functions:

- Over temperature alarm.
- Independent limit temperature alarm system.
- Leakage current and over-voltage protector.

Technical Specifications:

Model	QC-54-DI
Capacity (L)	54 L
Temperature Control Method	RT100
Temperature Range	RT+5~65°C
Temperature Stability	±1°C (at 37°C)
Temperature Uniformity	±2°C (at 37°C)
Ambient Temperature	5~30°C, recommend 25°C±2°C
Timing Range	0~999 h
Shelves	2pcs/adjustable
Power Consumption	345 W
Power Supply	AC220V±10%,50/60Hz
Internal Size (mm)	395*34*395 mm
External Size (W*D*H) mm	D type: 536*548*681 mm
Packing Size (W*D*H) mm	D type: 586*598*706 mm
Gross Weight (kg)	D type: 69 kg

Standard Accessories:

- Power Cord
- Standard Accessories
- User Manual