

QC – Low Speed Centrifuge Machine



FEATURES:

- Various rotor.
- Stepless speed regulation.
- Safety lid-lock, open lid automatically.
- Max. rotational speed: 500-5000 rpm.
- LCD display showing: RPM and RCF can be showed clearly.
- Large capacity: 8 or 6 of 15ml conical tubes with screw caps can be placed.

TECHNICAL DATA:

Model	QC-LC6S	QC-LC8S	QC-LC12S
Max speed	500 - 5000 r/min		500-4000 r/min
Max RCF	3074xg		2220xg

Max capacity	6*15ml (centrifuge tube)	8*15ml (centrifuge tube)	12*5/10ml (blood tube)
Speed regulation accuracy	± 30 r/min		
Timing range	10s ~ 99min59s		
Noise	≤ 60dB(A)		
Voltage	AC220V / 110V - 50Hz 5A		
Power	130W		
Dimensions	375*320*240 mm(L*W*H)		
Weight	5.5kg		
Rotor Pictures			