

DIGITAL ORBITAL SHAKER



Main Features

- 1.Speed: 40-200rpm for all models
- 2.Digital panel with speed and time control
- 3.Working pan size 25X25cm/32X32cm, with max load 5KG
- 4.Continuous operation and timed operation range 1min-23h 59min
- 5.DC brushless motor, safe, long life and free of maintenance. Operating silent and stable

Technical Specifications:

Model	QC-OS-105Pro
Motor Type	Brushless motor
Plate Size	250*250 mm
Turn	Circular
Operation Mode	20 mm
Max Load	5 kg
Speed Range	40 - 200 rpm
Timing Range	1min - 23h 59min
Display	LED display
Shell Paction Level	IP21
Power & Voltage	30W & 100-240V, 50/60Hz
Size	260*355*135mm
Certificate	ISO & CE Certificate