

Automatic Seed Counter Machine



Description:

Automatic seed Counter can automatically count many kinds of grains by setting number. In agriculture it is mainly used in grain count, exp. Rice, wheat, sorghum, corn, vegetables seed and so on.

Features:

- Two modes: Free count and bulk count;
- The counting speed is adjustable;
- High accuracy: Sensitivity levels for different size of particles;
- Automatic deceleration: when counting number is close to the setting number, the counting speed slow down automatically;
- Automatic motor Protection: the motor will stop work once it detects there is no particles on the plate for 60 seconds;
- User-friendly interface: LCD screen, mode selection, counting time, level and number of particles display clearly.
- Widely application: suitable for particles includes but is not limited to: crops seeds, vegetable seeds, feed, diamond, tablet.

Technical Specifications:

Model	QC-SC104
Count Range	1-12mm
Count Accuracy	±2‰
Count Speed	≥1000 /3min (Sesame)
Counting Range	1~99999
Bulk Counting	1-99999
Noise	<70db
Power	AC110V /60Hz and AC 220V /50Hz
Continuous working time	≥ 5Hr
Net Weight	5kg
Gross Weight	6.5kg
Packing Dimensions	400*350*260mm
Size of Machine	285*210*176mm
Operation Environment	-10°C- 50°C, <85% RH
Accessories	Standard accessories included.