


Adjustable Speed Vacuum Pump



Main Features

1. Working without pump oil.
2. Adjustable suction flow rate: 0-30 L/min.
3. Universal voltage 110V~220V.
4. Overload protection, over temperature protection, short circuit protection.
5. Can work without stopping for a long time, over 24 hours.

Technical Specifications:

| Model | QC-VP/30 |
|-----------------------------------------|--------------------------------|
| Type | Positive and Negative pressure |
| Pumping speed: (L/min) | 30 |
| Motor Speed (RPM) | 1500/Min. |
| Ultimate Vacuum/Negative Pressure (Mpa) | 0.085 |
| Max Pressure/Positive Pressure (bar) | 5.0 |
| Noise Level (DB) | <50 |
| Motor Power (W) | 50 |
| Diaphragm & Vales | NBR |
| Overload protection | Yes |
| Over temperature protection | Yes |
| Short circuit protection | Yes |
| Dimensions (mm) | 240x120x160 |
| Products weight (kgs) | 2.68 |
| Certificate | ISO & CE Certificate |