

# Auto Seed Blower Machine



## Description:

- To remove dead & light seeds and particles from the good seeds.
- Suitable for large seeds, including corn, wheat, rice, soybean, etc.

## Features:

- The shell is made of metal, strong and durable
- Lower noise, longer lifetime and power saving
- High-quality glass pipes, visual operation
- Automatic timing control
- Adjustable airflow according to particle size
- Better performance of the metal cover ventilation pipe

## Technical Specifications:

<b>Model</b>	<b>QC-BS1000</b>
Maximum wind pressure	680 pa
Maximum air volume	0.142 m <sup>3</sup> /min
Maximum wind speed	33 m/s
Speed	2800 rpm/min
Timing	0-9999 min
Cup mesh	28 mesh
Working noise	≤81dB
Motor power	550W
Working voltage	220V 50HZ
Dimensions	526mm×455mm×1680mm
Accessories	Standard accessories included.