

Digital Table Top Centrifuge Machine (QC-DC5T Plus)



Features:

- * All-metal structure multi-layer explosion-proof design, strong and durable.
- * Embedded microprocessor control, DC brushless motor drive, maintenance free, smooth operation.
- * 5-inch touch LCD screen, real-time display all running parameters, the operation interface is intuitive and simple.
- * Support small language customization,
- * Mechanical and electrical integration door lock, easy to use, safe and reliable.
- * Silicone integral sealing ring to avoid aerosol overflow.
- * With self-test, block, overspeed, imbalance, over temperature detection and alarm functions, excellent shock absorption device.
- * Store multiple groups of user-defined programs, boot for the last program.

Technical Specifications:

Max. Speed	5000 rpm	Noise	≤65dB (A)
Max. RCF	4807 × g	External Size (L*W*H)	510*630*420 mm
Max. Capacity	4*250ml	Net Weight	56 Kg
Speed Precision	±30 rpm	Gross Weight	70 Kg
Time Range	1min~99h59min	Supported Rotor	1.5ml, 2.0ml, 5ml, 10ml, 15ml, 50ml, 100 ml.
Power Supply	AC110/220V 50/60Hz	Certificate	ISO & CE Certified