

Sample Automatic Divider (QC-SD100)



Features:

- Extremely high division accuracy, Representative and reproducible division for accurate analysis results.
- The injection system and the sample separation system are integrated and easy to operate.
- Fully automatic design, preset time, speed and injection speed.
- The sampling chamber is optional, which is suitable for different sample.
- The split crown is made of aluminum material with abrasion resistance and difficult to adsorb samples.
- For samples with finer particle size, the cap is protected by dust cap.
- Split crown fastener structure design, easy disassembly.
- The overall design of the instrument is compact, convenient to clean.
- Low-noise drive.

Technical Specifications:

| | |
|-----------------------|------------------------------------|
| Sample Type | Bulk materials |
| Feed size | <12mm |
| Sample Divider batch | 5L |
| Sample Copies | 6/8/10 |
| Electrical Parameters | 220V/150W |
| Accessory | All Standard Accessories Included. |