

Vibrating Sieve Separator (QC-VSS50)



Features:

- Intermittent mode can be set to ensure the screening rate of the sample;
- With dry screen and wet screen two functions;
- No maintenance, very quiet operation;
- Unique three-dimensional throwing operation, to ensure the high precision of screening;
- LED5 inch touch screen, intelligent control, easy operation;
- Digital setting, can get accurate and repeatable results;
- Reliable screening results, comfortable and convenient humanized design;
- Can be connected to the computer to use special screening software for operation control, data processing record and analysis, screening evaluation.

Technical Specifications:

Size range of screening	20 microns to 25 mm	Storable parameter	20 groups
Amplitude range	0.2~3.0 mm	Maximum load	3.0 kg
Time frame	01 to 99 minutes	Rated power	400 W
Intermittent range	00~99 seconds	Rated voltage	220 V, 50 HZ
Drive screening mode	Electromagnetic drive	Instrument size	460*510*270 mm
Optional screen diameter	100\150\200\203 mm	Package size	440*670*680 mm
Screening number	Lv. 9	Instrument weight	48 kg