

Anaerobic Jar



Introduction

The Anaerobic Jar is a specialized laboratory container designed to create an oxygen-free environment for the cultivation of anaerobic microorganisms.

Features:

1. Excellent airtight performance, avoiding any undue leakage.
2. PMMA material, smooth surface, shock and impact resistant.
3. Jar and lid are both transparent, sealing between lid and jar consists of an excellent quality O-ring.
4. Strong jar-lid-clamp construction, guarantee the overall airtight performance.
5. Colored stainless-steel clamp with power electrostatic painting, durable and convenient for experimental management.
6. Adopting quick snap-shut coupling, enabling quick connection to the system in one swift move.
7. For BK-AJH series, configuration list: anaerobic jar, 15/25R rack, pressure gauge (standard: oil type, optional: oil-free type); For BK-AJG series, configuration list: anaerobic jar, 15/25 rack.

Technical Specification:

QC-AJL Series Pumping Type Jar: 3 to 5 minutes every time

Model	Contents	Petri Dish Holder	Package Size(W*D*H)	Gross Weight
QC-AJL15	1 stack of 6 dishes(ϕ 9~10cm)	15R	245*215*295mm	2.0kg
	ϕ =130mm, h=150mm; 1.5 liters			
QC - AJL25	1 stack of 12 dishes(ϕ 9~10cm)	25R	245*215*395mm	2.5kg
	ϕ =130mm, h=240mm; 2.5 liters			
QC - AJL35	1 stack of 18 dishes(ϕ 9cm)	15/25R	245*215*495mm	3.0kg
	ϕ =130mm, h=280mm; 3.5 liters			
QC -AJL50	3 stacks of 6 dishes(ϕ 9cm)	3/15R	370*370*340mm	5.6kg
	ϕ =240mm, h=120mm; 5 liters			
QC -AJL70	3 stacks of 8 dishes(ϕ 9cm)	3/15R	370*370*370mm	6.0kg
	ϕ =240mm, h=160mm; 7 liters			
QC -AJL100	3 stacks of 12 dishes(ϕ 9cm)	3/25R	370*370*470mm	7.0kg
	ϕ =240mm, h=230mm; 10 liters			

QC-AJI Series Hermetic Type Jar: About 4 hours every time

Model	Contents	Petri Dish Holder	Package Size(W*D*H)	Gross Weight
QC-AJI15	1 stack of 6 dishes(ϕ 9~10cm)	15R	245*215*295mm	2.0kg
QC-AJI25	1 stack of 12 dishes(ϕ 9~10cm)	25R	245*215*395mm	2.5kg
QC-AJI35	1 stack of 18 dishes(ϕ 9cm)	15/25R	245*215*495mm	3.0kg
QC-AJI50	3 stacks of 6 dishes(ϕ 9~10cm)	3/15R	370*370*340mm	5.6kg
QC-AJI70	3 stacks of 8 dishes(ϕ 9cm)	3/15R	370*370*370mm	6.0kg
QC-AJI100	3 stacks of 12 dishes(ϕ 9cm)	3/25R	370*370*470mm	7.0kg