

QC- Series Analytical Balance



Main Features

- Humanized design, novel design and five-side full-transparent windshield.
- A new-generation electromagnetic balance weighing sensor guarantees high precision.
- Application programs such as automatic fault detection, four-point linear calibration and over-load protection.
- Ultra-large LCD display brings clearer and more comfortable visual sense to operate.

Technical Specifications:

Model	FAQC-2204AB
Weighing Capacity (g)	0.0001 g - 220 g
Readability	0.0001 g (0.1 mg)
Repeatability	±0.0002 g
Linearity	±0.0003 g
Stabilization Time	Within 3 seconds
Pan Size	Φ80 mm
Gross Weight	22.1 lbs
Packaging Dimensions	13.57*18.08*17.49 inches
Sensor	Electromagnetic force sensor
Calibration Weight	200 g
Power	220V, 50Hz, AC
Warranty	12 Months
Type of calibration	ISO CAL
Certificate	ISO & CE Certificate

Standard Accessories:

1. Power Adapter
2. Weighting Standard
3. Weight Pan
4. User Manual