

## Digital Seed Purity Work-Board (DPW-1500 Plus)



### Features:

- The workbench is made of steel and wood, laminate surface
- Super-thin light panel, soft and even white light
- The polarized light against dazzling transmission observes the platform, electrically operated adjustable dark-field
- Transmission illumination function (The polarized light equipment can filter the background photo source completely the disturbance, may to the seed, the seedling, the plant leaf blade, the colony and so on careful observation examination, and is suitable in photographs the high-grade picture)
- Large diameter live arm type illumination magnifying glass.
- LED light

### Technical Specifications:

Magnification	0-7 X
Machine Size	130 * 65*75(cm)
Magnifier light power	20W
Observation table size	63 * 50(cm), power: 15W
Polarized light platform Size	27*27(cm), power; 16W
Accessory	Manual counter, grain-counting board