

QC – Glass Distilled Water Plant



FEATURES:

- **Economical Price.**
- **Glass construction with long life metal heating element.**
- **Two independent safety thermostats.**
- **Wall mountable.**
- **Electricity power saving.**
- **Easy to use.**
- **Duran Glassware system water distill unit.**

TECHNICAL DATA:

Model	QC-WDP/4L
Output	4 litres/hr. single distilled
pH	5.0 – 6.5
Conductivity, μScm^{-1}	3.0 – 4.0
Resistivity, mOhm-cm	0.25 – 0.3
Temperature	25 – 35°C
Pyrogen content	Pyrogen free
Water supply	1 litre/min 3 – 100 psi (200 – 700 kPa)
Electricity supply	220 or 240 V, 50-60 Hz, single phase
Power	3kW
Dimensions (WxDxH)	500 x 150 x 450 mm
IP Rating	31
Weight	8 kg
Warranty	1 Year

Note: Care is required to produce pyrogen free water and the output should be tested before use.