

QC- A Series Ducted Fume Hood



Main Features

- LED Display.
- Waterproof Socket.
- Water & Gas Remote Control.
- Double Layer (1.2mm sheet metal surface, melamine board inward).
- UV lamp for sterilization.
- Resistant to weak acid and alkali.
- Adjustable air speed: 9 levels
- Motorized front window, height adjustable.

- Microprocessor control system, LED display.
- With memory function in case of power-failure.

Technical Specifications:

Model	QC-FH1000(A)	QC-FH1200(A)	QC-FH1500(A)	QC-FH1800(A)
Internal Size (W*D*H)	830*680*740mm	1030*680*740mm	1330*680*740mm	1630*680*740mm
External Size (W*D*H)	1030*810*2200mm	1230*810*2200mm	1530*810*2200mm	1830*810*2200mm
Work Surface Height	850mm			
Max Opening	520mm			
Air Velocity	0.3~0.8m/s			
Noise	≤68dB			
Exhaust Duct	Φ300mm			
	PVC, Standard length: 4 meters.			
LED Lamp	8W*1	12W*1	16W*1	16W*1
UV Lamp	Emission of 253.7 nanometers for most efficient decontamination			
Blower	Built-in centrifugal blower; Speed adjustable.			
Front Window	5mm toughened glass; Motorized; Height adjustable.			
Power Supply	AC220±10%, 50/60Hz; 110V±10%, 60Hz			
Consumption	400W	400W	500W	500W
Material	Exterior	Cold-rolled steel with anti-bacteria powder coating.		
	Interior	High grade melamine board with good acid and alkali resistance function		
	Work Table	Chemical resistant phenolic resin		
Stand Accessory	V lamp*2, LED lamp, Water tap, Gas tap, Water sink, Base cabinet			
	Total load of 2 waterproof sockets: 500W			
	4-meter PVC exhaust duct, pipe strap*2			
Optional Accessory	Active carbon filter			
Gross Weight	286kg	317kg	350kg	400kg
Dimension (mm)	1190*990*1690mm	1390*990*1690mm	1690*990*1690mm	1990*990*1690mm