

## QC-Series Hot Air Oven



### Main Features

- External material: cold-rolled steel with anti-bacteria powder coating.
- Stainless steel inner chamber, round angle structure, adjustable shelves.
- Double-layer glass observation window, adjustable airtight buckle lock.
- PID microprocessor intelligent temperature controller.
- Over-temperature protection.

### Technical Specifications:

Model	QC-30F	QC-42F	QC-45F	QC-65F	QC-71F	QC-125F
Capacity (Liters)	30	42	45	65	71	125
Temp. Range	RT+10~300°C	RT+10~200/300°C	RT+10~300°C	RT+10~300°C	RT+10~300°C	RT+10~300°C
Temp. Precision	0.1°C					

Temp. Fluctuation	±1°C					
Ambient Temp.	5~40°C					
Timing Range	1~9999 min					
Shelves No.	2 pcs					
Display	LED					
Consumption Power	800W	850/1600W	1200W	1600W	1550/2400W	2300W
Power Supply	Standard: AC220V±10%, 50/60Hz					
Internal Size (W*D*H) mm	310*310*310	350*350*350	350*350*350	400*360*450	450*352*450	500*450*550
External Size (W*D*H) mm	460*510*695	630*630*552	500*550*735	550*550*840	735*630*652	636*680*915
Gross Weight (kg)	34	60	39	47	85	58