

QC- Series Infrared Sterilizer



Main Features

- Sterilization only takes 5-7 seconds
- Small volume, light weight
- Easy to clean
- Used in anaerobic chamber
- Wide range of applications (inoculating loop, inoculating needle or culture tube mouth)

Technical Specifications:

Model	QC-S22-A	QC-S22-B
Power	220V/150W	220V/240W
The Center High-Temperature	825°C ± 50°C	
Standby Constant Temperature	Heating with two gears 400°C / 850°C	
The Max-Diameter of Sterilizer Goods	φ 14mm	φ 35mm
The Length of Sterilizer	150mm	100mm
Heating Time	10min	
Working Condition Temperature	5 ~ 40°C	
Relative Humidity	≤90%	
Dimensions	175*130*185mm	147*130*200mm
Net Weight	1.3kg	1.4kg