

## Seed Purity Work-Board (SPW-1300 Plus)



### Features:

- Both sides of the work table adopts slope design according to the human body engineering principle, the operation is comfortable and labor-saving.
- Removable drawer is under the instrument to store well-chosen seeds.
- Extended edge is designed on the left, right and back side of workbench to avoid seeds' lateral leakage.
- Glass surface can be drawled out to replace tube easily.
- The original seeds-fall holes are advantageous to take out the chosen seeds.
- Metal structure, ultra-thin observation platform, provided with desktop illumination magnifying glass.

### Technical Specifications:

Magnification	0-7 times
Power	≤30W, 220V/50Hz
Surface size	320 * 320mm
Table base size	930 * 330mm
Luminous panel size	260*260mm