

## QC Series-Sieve Shaker



### Main Features

- Screen mesh size range: from 4 mesh to 500 mesh.
- Timing range: 0~99min59s, time can be set from seconds to hours freely.
- Typically requires 50~150g of sample material.

### Technical Specifications:

|                       |  |
|-----------------------|--|
| Model                 | QC-SS 400                              |
| Working Layer         | 1~8 layers                             |
| Screen Mesh Size      | 10, 20,50, 80,100 mesh etc.as required |
| Diameter              | 200mm                                  |
| Material              | Stainless steel 304                    |
| Timing Range          | 0~99min 59 sec                         |
| Speed                 | 1400r/min                              |
| Power Consumption     | 120W                                   |
| Power Supply          | AC110/220V±10%, 50/60Hz                |
| External Size (W*D*H) | 400*330*370mm                          |
| Package Size (W*D*H)  | 540*490*460mm                          |
| Gross Weight          | 47kg                                   |