

Digital Ultrasonic Cleaner Bath



Main Features

- Structural design: The shell of the whole ultrasonic cleaner is formed by laser cutting of SUS304 stainless steel, and the surface adopts anti-fingerprint and oleophobic treatment.
- Display and control design: 4.3-inch backlit LCD display, clear and bright in three colors; sensitive and smooth response for touch control operation.
- Mute design: UVF intelligent mute series is born in the name of mute, and the best multiple sound reduction design truly ensures that it is 26 decibels lower than similar products.
- Circuit design: Built-in German IBDI ultrasonic generator chip can realize ultrasonic circuit sound energy conversion, with frequency sweep function, high frequency stability, and industrial-grade panels.
- Appearance design: Good-looking glass panel design, the whole machine is colorful and bright, and has a scientific and technological texture.

Time Control (Min)	1-999	1-999	1-999	1-999	1-999	1-999	1-999	1-999	1-999	1-999
Internal tank size (mm)	240x1	300x	300x	300x	300x	330x	500x	500x	500x	500x
	40x	150x	240x	240x	240x	300x	300x	300x	300x	300x
	100	150	150	150	150	150	150	150	150	200
External tank size (mm)	264x	325x	328x	328x	328x	355x	530x	530x	530x	530x
	162x	175x	266x	266x	266x	326x	328x	328x	328x	328x
	242	292	290	290	290	292	292	292	292	342
Standard Configuration	Sound reduction cover, silicone rubber gasket, Bracket, Drain Set									