

QC – Veterinary Urine Analyzer



FEATURES:

- Program storage is 64M+512K, Movement memory is 32M+96K.
- Single chip control the whole system which enhance the stability and reliability.
- Color date to identify, not gray or dyeing recognition. Avoid the dark color urine influence (interference).
- Printer inside, small and big formats printing styles.
- Connect with urine sediment analyzer and other instruments.
- Test results realize the directly comparison.
- Barcode reader system, scan and store the test results.
- Any animal urine can be tested.
- Testing cost very low.

TECHNICAL DATA:

Model	QC-VET UA-01
Display	5-6-inch (touching) color LCD screen
Test parameters	URS-10T, URS-11, URS-12
wavelength	558nm, 635nm, 720nm
Test item	Leu, Nit, Uro, Pro, pH, Blood, SG, VC, Ket, Bil, Glu, Microalbumin
Test speed	Up to 120 tests/hour at Continuous
Memory Size	1000 test results
Option	10 or 11 parameters
Print	Thermal printer inside, with outer stylus printer interface
Host-communication	PC connected with RS232
Language	English, European, Korean, Japanese, China
Environment Condition	0-40°C. RH< 85%
Power Supply	AC 220V (=15%), 50-60Hz
Power	60W
Dimensions (W/D/H)	350*320*150mm
Weight	9 Kg